



Audiological Service & Supply Company

800-755-8272

Audiometry.com

Optec 6500

**Patented Controlled Lighting
for the most accurate color
vision testing available**



Features:

- New Unsurpassed Homogenous Illumination
- All Light levels continually calibrated and microprocessor controlled
- Binocular Glare Testing
- Target illumination with 3.0 cd/m for night or 85 cd/m for day with a value range of 1.0 cd/n to 125 cd/m
- Glare luminance at distance with 1 Lux for night glare, and 10 Lux for day glare with a glare luminance value range of 1 Lux to 50 lux
- Contrast Sensitivity test sequence to evaluate the patient's overall functional vision. Choose from either Sine Wave (F.A.C.T.) or Letter (C.A.T.)
- Potential Acuity results with a quick assessment of macular function in cataract patients that documents if surgery is expected to improve visual performance.

800-755-8272



Fax:773-4274863